

# Wachsumrechnungstabelle

CE 0197

Artikel-Nr.	Dentallegierung	Wachsgewicht in Gramm inklusive Gusskanalsystem (ohne Gusskegel)																										
		0,2	0,3	0,4	0,5	0,6	0,7	0,8	0,9	1,0	1,2	1,4	1,6	1,8	2,0	2,2	2,4	2,6	2,8	3,0	3,2	3,4	3,6	3,8	4,0	4,2	4,4	4,6

Bio-Aufbrennlegierungen		Verblendbar mit Noritake EX-3																												
1001	Bio Keramikgold G	4	6	8	10	12	14	15	17	19	23	27	31	34	38	42	46	49	53	57	61	65	68	72	76	80	83	87	91	95
1002	Bio Keramikgold TI	4	6	8	10	11	13	15	17	19	23	26	30	34	38	41	45	49	52	56	60	64	67	71	75	78	82	86	90	93

Hochgoldhaltige Aufbrennlegierungen		Verblendbar mit Noritake EX-3																												
1003	Keramikgold U	4	5	7	8	10	12	13	15	17	20	23	26	29	33	36	39	42	46	49	52	55	59	62	65	68	72	75	78	81
1004	Keramikgold P6	4	6	8	10	11	13	15	17	19	23	26	30	34	38	41	45	49	52	56	60	64	67	71	75	78	82	86	90	93
1008	Keramikgold H	4	6	8	10	11	13	15	17	19	23	26	30	34	38	41	45	49	52	56	60	64	67	71	75	78	82	86	90	93
1005	Europe 76	4	6	8	10	12	14	15	17	19	23	27	31	34	38	42	46	49	53	57	61	65	68	72	76	80	83	87	91	95

Goldreduzierte Aufbrennlegierungen		Verblendbar mit Noritake EX-3																												
1101	Palladiumgold 52	3	5	6	7	9	10	12	13	15	17	20	23	26	29	32	34	37	40	43	46	49	51	54	57	60	63	66	68	71
1102	Palladiumgold 6	3	4	5	6	7	8	10	11	12	14	17	19	21	24	26	28	30	33	35	37	40	42	44	47	49	51	54	56	58
1103	Palladiumgold S1	3	4	5	6	7	8	9	10	12	14	16	18	21	23	25	27	30	32	34	36	39	41	43	46	48	50	52	55	57
1104	Europe 45	3	4	5	6	7	8	9	10	12	14	16	18	21	23	25	27	30	32	34	36	39	41	43	46	48	50	52	55	57

Universallegierungen (Guss- oder Aufbrennlegierungen)		Verblendbar mit Quattro Ceramic HI <sup>®</sup>																												
1201	Universalgold Norm	4	5	7	9	10	12	14	15	17	20	24	27	30	34	37	40	44	47	50	54	57	60	64	67	70	74	77	80	84
1202	Universalgold PDF	4	5	7	9	10	12	14	15	17	20	24	27	30	34	37	40	44	47	50	54	57	60	64	67	70	74	77	80	84
1207	Universalgold PDF Pur	4	5	7	9	10	12	14	15	17	20	24	27	30	34	37	40	44	47	50	54	57	60	64	67	70	74	77	80	84
1203	Universalgold KF	3	5	7	8	10	11	13	14	16	19	22	25	29	32	35	38	41	44	47	51	54	57	60	63	66	69	73	76	79
1204	Universalgold RE	3	4	6	7	9	10	11	13	14	17	19	22	25	28	30	33	36	39	41	44	47	50	52	55	58	61	63	66	69
1206	Universalgold Eco	3	4	5	6	8	9	10	11	13	15	18	20	23	25	28	30	33	35	38	40	43	45	48	50	53	55	58	61	63
1205	Universalgold P	2	3	4	5	6	8	9	10	11	13	15	17	19	21	23	25	27	29	31	33	35	37	39	42	44	46	48	50	52

Goldgusslegierungen																														
2005	Inlaygold T	3	5	6	8	10	11	13	14	16	19	22	25	28	31	34	37	40	43	47	50	53	56	59	62	65	68	71	74	77
2001	Biogussgold	4	6	8	9	11	13	15	17	18	22	26	29	33	37	40	44	47	51	55	58	62	65	69	73	76	80	84	87	91
2002	Gussgold M	3	5	6	8	10	11	13	14	16	19	22	25	28	31	34	37	40	43	47	50	53	56	59	62	65	68	71	74	77
2003	Gussgold P 4	3	5	7	8	10	11	13	14	16	19	22	25	28	32	35	38	41	44	47	50	53	56	60	63	66	69	72	75	78
2004	Gussgold KF	3	5	7	8	10	11	13	14	16	19	22	25	28	32	35	38	41	44	47	50	53	56	60	63	66	69	72	75	78
2006	Hybridgold	4	6	8	10	12	14	16	18	20	23	27	31	35	39	43	47	50	54	58	62	66	70	74	78	81	85	89	93	97
2101	Gussgold RE	3	4	6	7	9	10	11	13	14	17	20	22	25	28	31	33	36	39	42	44	47	50	53	56	58	61	64	67	69
2102	Gussgold RCF	3	4	6	7	9	10	11	13	14	17	20	22	25	28	31	33	36	39	42	44	47	50	53	56	58	61	64	67	69

Nettolegierungen		Verblendbar (siehe WAK-Wert)																												
1105	Palladium N	3	4	5	6	7	8	9	10	12	14	16	18	21	23	25	27	30	32	34	36	39	41	43	46	48	50	52	55	57
1209	Unigold N	2	3	4	5	6	8	9	10	11	13	15	17	19	21	23	25	27	29	31	33	35	37	39	42	44	46	48	50	52
1210	Kombigold N	3	4	6	7	9	10	11	13	14	17	20	22	25	28	31	33	36	39	42	44	47	50	53	56	58	61	64	67	69



GOLDQUADRAT GmbH  
 Büttnerstraße 13 · 30165 Hannover  
 Tel.: +49 (0)511 449897-0 · Fax: 449897-44  
 info@goldquadrat.de · www.goldquadrat.de



Stand: April 2010

